

Product datasheet for **TP327656L**

DHRS1 (NM_001136050) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dehydrogenase/reductase (SDR family) member 1 (DHRS1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227656 protein sequence Red =Cloning site Green =Tags(s)

MAAPMNGQVCVVTGASRGIGRGIALQLCKAGATVYITGRHLDTLRVAQEAQSLGGQCVWCDSSQESE
VRSLFEQVDREQQGRDLVNNAYAGVQTLNTRNKAFWETPASMWDDINNVGLRGHYFCSVYGARLMVP
AGQGLIVISSPGSLQYMFNVLYGVGKAACDKLAADCAHELRRHGVSCVSLWPGIVQTELLKEHMAKEEV
LQDPVLKQFKSAFSSAETTELSGKCVVALATDPNILSLSGKVLPSCDLARRYGLRDVDGRPVQDYLSLSS
VLSHVSGLGWLASYLPSFLRVPKWIILYTSKF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	33.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129522
Locus ID:	115817



[View online »](#)

UniProt ID: [Q96LJ7](#)

RefSeq Size: 1495

Cytogenetics: 14q12


RefSeq ORF: 939

Synonyms: SDR19C1

Summary: This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Nov 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified DHRS1 protein (Cat# [TP327656]). The protein was produced from HEK293T cells transfected with DHRS1 cDNA clone (Cat# [RC227656]) using MegaTran 2.0 (Cat# [TT210002]).