

Product datasheet for TP327659L

NECAP2 (NM_001145277) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NECAP endocytosis associated 2 (NECAP2), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227659 representing NM_001145277 Red =Cloning site Green =Tags(s)

MEESGYESVLCVKPDVHVYRIPPRATNRGYRAAEWQLDQPSWSGRLRITAKGQMAYIKLEDRTSGELFAQ
APVDQFPGTAVESVTDSSRYFVIRIEDGNGRRRAFIGIFGDRGDAFDNFVALQDHFVWVKQQCEFAKQAQ
NPDQGPKLDLGFKEGQTIKLNIANMKKKEGAAGNPRVRPASTGGLSLLPPPPGGKTSTLIPPPGEQLAVG
GSLVQPAVAPSSDQLPARPSQAQAGSSSDLSTVFPVHTSGKALPHLGQRKEDEALLSWPVFGA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	29.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138749
Locus ID:	55707



[View online »](#)


UniProt ID: [Q9NVZ3](#)

Cytogenetics: 1p36.13

RefSeq ORF: 819

Summary: This gene likely encodes a member of the adaptin-ear-binding coat-associated protein family. Studies of a similar protein in rat suggest a role in clathrin-mediated endocytosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]

Product images:



Coomassie blue staining of purified NECAP2 protein (Cat# [TP327659]). The protein was produced from HEK293T cells transfected with NECAP2 cDNA clone (Cat# [RC227659]) using MegaTran 2.0 (Cat# [TT210002]).