

Product datasheet for **TP327670**

BSDC1 (NM_001143890) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BSD domain containing 1 (BSDC1), transcript variant 4, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227670 representing NM_001143890 Red =Cloning site Green =Tags(s)
	 MAEGEDVGWWRSLQQSYQAVKEKARLYSLQSDPATYCNEPDGPPPELFDWLSQFCLEEKKGISELLVG SPSIRALYTKMVPAAVSHSEFWHRYFYKVHQLQEQARRDALKQRAEQSISEEPGWEEEEELMGISPIS PKEAKVPVAKISTFPEGEPGPQSPCEENLVTSVEPPAEVTPSESSISLVTQIANPATAPEARVLPKDL SQKLEASLEEQGLAVDVGETGSPPIHSKPLTPAGHTGGPEPRPPARVETLREEAPTDLRVFELNSDSG KSTPSNNGKKGSSTDISDWEKDFDLDMTEEEVQMALSKVDASGELEDVEWEDWE TRTRPLEQLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	37 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001137362</u>
Locus ID:	55108



[View online »](#)

UniProt ID: [Q9NW68](#)
Cytogenetics: 1p35.1
RefSeq ORF: 1005

Product images:

Coomassie blue staining of purified BSDC1 protein (Cat# TP327670). The protein was produced from HEK293T cells transfected with BSDC1 cDNA clone (Cat# [RC227670]) using MegaTran 2.0 (Cat# [TT210002]).