

Product datasheet for **TP327688M**

TRAPPC3L (NM_001139444) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BET3 like (<i>S. cerevisiae</i>) (BET3L), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227688 representing NM_001139444 Red =Cloning site Green =Tags(s)
	MSRPAHRRPEYHKINKDLFVLTYGALVAQLCKDYEKDEDVNQYLDKMGYGIGTRLVEDFLARSCVGRCHS YSEIIDIIAQVAFKMYLGITPSVTCNNSSKNEFSLILEKNPLVEFVEELPAGRSSLCYCNLLCGIIRGAL EMVHLAADVTFLQDRLKGDSVTEIGITFLKKRDEKKYRGKK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001132916</u>
Locus ID:	100128327
UniProt ID:	<u>Q5T215</u>
Cytogenetics:	6q22.1



[View online »](#)

RefSeq ORF:	543
Synonyms:	bA259P20.2; BET3L
Summary:	May play a role in vesicular transport from endoplasmic reticulum to Golgi.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TRAPPC3L protein (Cat# [TP327688]). The protein was produced from HEK293T cells transfected with TRAPPC3L cDNA clone (Cat# [RC227688]) using MegaTran 2.0 (Cat# [TT210002]).