

Product datasheet for **TP327752M**

SDK2 (NM_001144952) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sidekick homolog 2 (chicken) (SDK2), 100 µg
Species:	Human
Expression Host:	HEK293T



[View online »](#)

Expression cDNA Clone or AA Sequence: >RC227752 representing NM_001144952
Red=Cloning site Green=Tags(s)


MWGLLIWTLALHQAARAQDDVSPYFKTEPVRTQVHLEGNRLVLTCAEGSWPLEFKWLHNNRELTK
 F
 SLEYRYMITSLDRTHAGFYRCIVRNRMGALLQRQTEVQVAYMGSFEEGEKHQSVSHGEAAVIRAPRIASF
 PQQPQVTWFRDGRKIPSSRIAITLENTLVILSTVAPDAGRYVQAVNDKNGDNKTSQPITLTVENVGGA
 DPIAPTIIPKNTSVVAGTSEVTECVANARPLIKLHIIWKKDGVLLSGGISDHNRLTIPNPTGSDAG
 YYECEAVLRSSSVPSVVRGAYLSVLEPPQFVKEPERHITAEMEKVVDIPCQAKGVPPPSITWYKDAAVE
 VEKLTRFRQRNDGGLQISGLVDDTGMFQCFARNAAGEVQTSTYLAVTSIAPNITRGPLDSTVIDGMSV
 LACETSGAPRPATWQKGERILASGSVQLPRFTPLESGSLISPTHSIDAGTYTCLATNSRGVDEASADL
 VVWARTRITKPPQDQSVIKGTQASMVCGVTHDPRVTIRYIWEKDGTATLGESHPRILDRNGSLHISQTW
 SGDITGTYTCRISAGGNDRSASHLRVRQLPHAPEHPVATLSTVERRAINLTWTKPFDGNSPLIRYILEMS
 ENNAPWTVLLASVDPKATSVTVKGLVPARSYQFRLCAVNDVGKQGFSKDTERVSLPEEPPTAPPQNVIAS
 GRTNQSIMIQWQPPPEHQNGILKGYIIRYCLAGLPVGYQFKNITDADVNNLLLEDLIWNTNIEIEVAAY
 NSAGLGVSYSKVTEWTLQGVPTVPPGNVHAEATNSTTIRFTWNPSPQFINGINQGYKLIawePEQEEV
 TMVTARPNFQDSIHVGFVSGLKKEFTEYFTSVLCFTTPGDGPRSTPQLVRTHEDVPGVPVGHLSFSEILDTS
 LKVSWEQEPGEKNGILTYRISWEEYNRTNTRVTHYLPNVTLEYRVTGLTALTTYIEVAAMTSKGGQGV
 ASTISSGVPELPGPPTNLGISNIGPRSVTLQFRPGYDGKTSISRWLVEAQGVVGEGEWLLIHQLSNE
 PDARSMVDPDLNPFTCYSFRMRQVNIVGTSPPSQPSRKIQLTQAPPDMPANVSLRTASETSLWLRWMP
 L
 PEMEYNGNPESVGYKIKYSRSDGHGKTLSHVQDRVERDYTIEDLEEWTEYRVQVQAFNAIGSGPWSQTV
 VGRTRVSSSGPTNVSALETTSSSMLVRWSEVPEADRNGLVLGYKVMYKEKSDTQPRFWLVEGNSSRS
 AQLTGLGKYVLYEVQVLAFTRIGDGSPPHPIERTLDDVPGPPMGILFPEVRTTSVRLIWQPPAAPNGI
 ILAYQITHRLNTTANTATVEVLAPSARQYTATGLKPESVYLFRITAQTRKGWGEAAEALVTTTEKDRP
 QPPSRPMVQQEDVRARSVLLSWEPSDGLSPVRYTYITRELPSGRWALHSASVSHNASSFIVDRPKFT
 SYKFRVKATNDIGDSEFSESESLTTLQAAPDEAPTILSVTPHTTTSVLIRWQPPAEDKINGILLGFRIR
 YRELLYEGLRGFTLRGINNPGATWAEELTSMYSMRNLSRPSLTQYELDNLNKHRRYEIRMSVYNAVGEGPS
 SPPQEVFVGEAVPTAAPRNVVHVGATATQLDVTWEPPPLDSQNGDIQGYKIYFWEAQRGNLTERVKTFL
 AENSVKLKNTGYTAYMVSVAAFNAAGDGRSTPTQGGTQQAAPSAPSSVKFSELTTSVNVSWEPQFP
 NGILEGYRLVYEPSPVDGVSKIIVTDVKGNLPLWLVKDLAEGVTYRFRIRAKTFTYGPEIEANVTGP
 GEGAPGPPGVPIIVRYSSAIAIHWSGDPGKGPITRYVIEARPSDEGLWDILIKDIPKEVSSYTFMSDIL
 KPGVSYDFRVIANDYGFPTSSPSQSVPAQKANPFYEEWWFLVIALVGLIFILLVFLIIRGQSKKY
 AKKTDSGNSAKSGALGHSEMMSLDESSFALELNRRLSVKNSFCRKNGLYTRSPRPSGSLHYSDDEV
 TKYNDLIPAESSSLTEKPSEISDSQGSSEYEVDSNHQKAHSFVNHYISDPTYNSWRRQKQGISRAQAY
 SYTESDSGEPDHTTVTNSTSTQQGSLFRPKASRTPTPQNPPNPPSQQSTLYRPPSSLAPGSRAPAGFSS
 FV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 239.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001138424</u>
Locus ID:	54549
UniProt ID:	<u>Q58EX2</u>
Cytogenetics:	17q25.1
RefSeq ORF:	6516
Summary:	The protein encoded by this gene is a member of the immunoglobulin superfamily. The protein contains two immunoglobulin domains and thirteen fibronectin type III domains. Fibronectin type III domains are present in both extracellular and intracellular proteins and tandem repeats are known to contain binding sites for DNA, heparin and the cell surface. This protein, and a homologous mouse sequence, are very similar to the Drosophila sidekick gene product but the specific function of this superfamily member is not yet known. Evidence for alternative splicing at this gene locus has been observed but the full-length nature of additional variants has not yet been determined. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SDK2 protein (Cat# [TP327752]). The protein was produced from HEK293T cells transfected with SDK2 cDNA clone (Cat# [RC227752]) using MegaTran 2.0 (Cat# [TT210002]).