

## Product datasheet for TP327796

### CCDC169 (NM\_001144981) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 13 open reading frame 38 (C13orf38), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227796 representing NM_001144981 Red=Cloning site Green=Tags(s)

MKEERNYNFDGVSTNRLKQQLLEEVRRKKDAVQLSIFELRHKITELEAKLNTDNEGSEWKTRYETQLELND  
ELEKQIVYLKEKVEKIHGNSSDRLSSIRVYERMPVESLNTLLKQLEEEKKTLESQVKYYALKLEQESKAY  
QKINNERRTYLAEMSQGSGLHQVSKRQQVDQLPRMQENLVKTGRYNPAKQKTVSAKRGVPVKITRPNHLP  
ELHP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	25.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001138453</a>
Locus ID:	728591



[View online »](#)

UniProt ID: [A6NNP5](#)  
Cytogenetics: 13q13.3  
RefSeq ORF: 642  
Synonyms: C13orf38

### Product images:



Coomassie blue staining of purified CCDC169 protein (Cat# TP327796). The protein was produced from HEK293T cells transfected with CCDC169 cDNA clone (Cat# [RC227796]) using MegaTran 2.0 (Cat# [TT210002]).