

## Product datasheet for **TP327808**

### Transcription Factor SP9 (SP9) (NM\_001145250) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens Sp9 transcription factor homolog (mouse) (SP9), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227808 representing NM_001145250 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MATSILGEEPRFGTTPLAMLAATCNKIGNTSPLTTLPESSAFAKGGFHPWKRSSSSCNLGSSLSGFAVAT GGRGSGGLAGGSGAANSFAFLASTSPTSSAFSSDYGGLFSNSAAAAAAGVSPQEAGGQSAFISKVHTT AADGLYPRVGMHPYESWYKSGFHSTLAAGEVTNGAASSWWDVHSSPGSWLEVQNPAGGLQSSLHSGAPQ ASLHSQLGTYNPDFSSLTHSAFSSTGLGSSAAAASHLLSTSQHLLAQDGFKPVLPYSYDSSAAVAAAAAAS AMISGAAAAAAGGSSARSARRYSGRATCDPCNCQEAERLGPAGASLRRKGLHSCHIPGCGKVYKGTSHLK AHLRWHTGERPFVCNWLFCGKRFRTRSDQLRHLRTHTEKRFACPVCKRFRMSDHLKHIKTHNGGGGG KKGSDSDTDASNLETPRSESPDLILHDSGVSAARAAAAAAAAAAAAAAAAASAGGKEAASGPND  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	48.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_001138722](#)

Locus ID: 100131390

UniProt ID: [P0CG40](#)


Cytogenetics: 2q31.1

RefSeq ORF: 1452

Synonyms: ZNF990

**Summary:** Transcription factor which plays a key role in limb development. Positively regulates FGF8 expression in the apical ectodermal ridge (AER) and contributes to limb outgrowth in embryos (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified SP9 protein (Cat# TP327808). The protein was produced from HEK293T cells transfected with SP9 cDNA clone (Cat# [RC227808]) using MegaTran 2.0 (Cat# [TT210002]).