

Product datasheet for **TP327854**

CD133 (PROM1) (NM_001145847) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human prominin 1 (PROM1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227854 protein sequence Red =Cloning site Green =Tags(s)

MALVLGSLLLLGLCGNSFSGGQPSSTDAPKAWNYELPATNYETQDSHKAGPIGILFELVHIFLYVWQPRD
 FPEDTLRKFLQKAYESKIDYDKIVYYEAGIILCCVLGLLFIILMPLVGYFFCMCRCCNKCGGEMHQRQKE
 NGPFLRKCF AISLLVICIIISIGIFYGFVANHQVRTRIKRSRKLADSNFKDLRLLNETPEQIKYILAQY
 NTKDKAFTDLNSINSVLGGGILDRLRPNIIPVLDEIKSMATAIKETKEALENMNSTLKSLHQQSTQLSS
 SLTSVKTSLRSSLNDPLCLVHPSETCNSIRLSLSQLNSNPELRQLPPVDAELDNVNNVLRDLGLVQQ
 GYQSLNDIPDRVQRQTTTAVAGIKRVLNSIGSDIDNVTQRLPIQDILSAFSVYVNNNTESIHRNLPTLEE
 YDSYWWLGGVLICSLTLVIFYLGLLCGVCYDRHATPTTRGCVSNTGGVFLMVGVLGSLFCWILMI
 IVLTFVFGANVEKLICEPYTSKELFRVLDTPYLLNEDWEYLSGKLFNKSKMMLTFEQVYSDCKKNRGT
 YGTLHLQNSFNISEHLNINEHTGSISSELESKLVNLNIFLLGAAGRKNLQDFAACGIDRMNYDSYLAQTG
 KSPAGVNLFSFAYDLEAKANSLPPGNLRNSLKRDAQTIKTIHQQRVLPQSLSTLYQSVKILQRTGNGL
 LERVTRILASLDFAQNFITNNTSSVIIIEETKKYGRITIIIGYFEHYLQWIEFSISEKVASCKPVATALDTAV
 DVFLCSYIIDPLNLFWFGIGKATVFLLPALIFAVKLAKYYRRMDESDVYDDVETIPMKNMENGNGYHKD
 HVYGIHNPVMTSPSQH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	96.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001139319
Locus ID:	8842
UniProt ID:	O43490 , A0A0A0N0M1
RefSeq Size:	4257
Cytogenetics:	4p15.32
RefSeq ORF:	2568
Synonyms:	AC133; CD133; CORD12; MCDR2; MSTP061; PROML1; RP41; STGD4
Summary:	This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation. Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Product images:



Coomassie blue staining of purified PROM1 protein (Cat# TP327854). The protein was produced from HEK293T cells transfected with PROM1 cDNA clone (Cat# [RC227854]) using MegaTran 2.0 (Cat# [TT210002]).