

Product datasheet for TP328058

LOC728392 (NM_001162371) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC728392 (LOC728392), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228058 representing NM_001162371 Red =Cloning site Green =Tags(s) MADGLFRRRPWGLEQIRPDPESEGLFDKPPPEDPPAARGPRASAAAGKKAGRRAGGRAQGGGRAGQPPKA A SRPPPKKEAPPLDEGCYLDHFPHLISIFIYAAIAFSITSCIFTYIHLQLA TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	12.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001155843
Locus ID:	728392
UniProt ID:	A0A494C111
Cytogenetics:	17p13.2


[View online »](#)

RefSeq ORF: 357

Product images:



Coomassie blue staining of purified LOC728392 protein (Cat# TP328058). The protein was produced from HEK293T cells transfected with LOC728392 cDNA clone (Cat# [RC228058]) using MegaTran 2.0 (Cat# [TT210002]).