

Product datasheet for **TP328152**

ANKRD66 (NM_001162435) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ankyrin repeat domain containing gene (LOC100287718), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC228152 representing NM_001162435
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MHSPGGLTTTKMAGPLPRVSDSLFSAMELAKMSDMTKLHQAVAAGDYSLVKKILKKGCLDPNYKDVDWND
 RTPLHWAAIKGQMEVIRLLIEYGARPCLVTSVGWTPAHFAAEAGHLNLIKTLHALHAAIDAPDFFGDTPK
 RIAQIYGQKACVAFLEKAEPECQDHRCAAQKGLPLDERDEDWDAKKRELESLPSLNQNMNKKKKS
 PTRPSNTKGRRV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	24.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001155907</u>
Locus ID:	100287718
UniProt ID:	<u>B4E2M5</u>




[View online »](#)

Cytogenetics: 6p12.3

RefSeq ORF: 666

Product images:



Coomassie blue staining of purified ANKRD66 protein (Cat# TP328152). The protein was produced from HEK293T cells transfected with ANKRD66 cDNA clone (Cat# [RC228152]) using MegaTran 2.0 (Cat# [TT210002]).