

## Product datasheet for TP328247M


### TSPY4 (NM\_001164471) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens testis specific protein, Y-linked 4 (TSPY4), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228247 representing NM_001164471 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MRPEGSLTYRVPRLRQGFVGRVAAQALVCASAKEGTAFRMEAVQEGAAGVESEQAALGEEAVLLDDI MAEEVVAEEVVAEEGLVERREEAQAQAVPGPGPMTPEALEELLAVQVELEPVNAQARKAFSRQR EKMERRRKPRLDRRGAVIQSVPGFWANVIANHPQMSALITDEDEDMLSYMVSLEVEEEKHPVHLCKIMLF FRSNPYFQNKVITKEYLVNITEYRASHSTPIEWYPDYVEAYRRRHNNSSLNFFNWFSDFHNFAGSNKIAE ILCKDLWRNPLQYYKRMKPPEEGTETSGDSQLLS  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	35.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001157943</a></u>
Locus ID:	728395


[View online »](#)

UniProt ID: P0CV99  
 Cytogenetics: Yp11.2  
 RefSeq ORF: 942  
 Synonyms: TSPY8; TSPY10

**Product images:**


Coomassie blue staining of purified TSPY4 protein (Cat# [TP328247]). The protein was produced from HEK293T cells transfected with TSPY4 cDNA clone (Cat# [RC228247]) using MegaTran 2.0 (Cat# [TT210002]).