

Product datasheet for TP328314

SKA3 (NM_001166017) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human spindle and kinetochore associated complex subunit 3 (SKA3), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228314 representing NM_001166017 Red =Cloning site Green =Tags(s)
	MDPIRSFCGKLRSLASTLDCETARLQRALDGEESDFEDYPMRILYDLHSEVQTLKDDVNILLDKARLENQ EGIDFIKATKVLMEKNSMDIMKIREYFQKYGYSPRVKKNSVHEQEAINSDELSNCENFQKTDVKDDLSD PPVASSCISEKSPRSPQLSDFGLERYIVSQVLPNPPQAVNNYKEEPVIVTPPTKQSLVKVLKTPKCALKM DDFECVTPKLEHFGISEYTMCLNEDYTMGLKNARNNKSEEAIDTESRLNDNVFATPSPIIQQLEKSDAEY TNSPLVPTFCTPGLKIPSTKNSIALVSTNYPLSKTNSSSNDLEVEDRTSLVLNSDTCFENLTDPSPTIS SYENLLRTPTPPEVTKIPEDILQKFQWIYPTQKLNKMR
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	43.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001159489</u>



[View online »](#)

Locus ID: 221150

UniProt ID: [Q8IX90](#)


Cytogenetics: 13q12.11

RefSeq ORF: 1164

Synonyms: C13orf3; RAMA1

Summary: This gene encodes a component of the spindle and kinetochore-associated protein complex that regulates microtubule attachment to the kinetochores during mitosis. The encoded protein localizes to the outer kinetechore and may be required for normal chromosome segregation and cell division. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]

Product images:



Coomassie blue staining of purified SKA3 protein (Cat# TP328314). The protein was produced from HEK293T cells transfected with SKA3 cDNA clone (Cat# [RC228314]) using MegaTran 2.0 (Cat# [TT210002]).