

Product datasheet for **TP328385L**

FAM90A5 (NM_001164455) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens family with sequence similarity 90, member A5 (FAM90A5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228385 representing NM_001164455 Red =Cloning site Green =Tags(s)

MMARRDPTSWAKRLVRAQTLQKQRRAPVGPRAPPPDEEDPRLKCKNCGAFGHTARSTRCPMKCWKAA
 LVP
 ATLGKKEGKENLKPWKPRVEANPGPLNKDKGEKEERPRQQDPQRKALLHMFSGKPPEKPLPNGKGSTEP
 S
 DYLRVASGPMVPVHTTSKRPRVDPVLADGSATEMSDRGSLASLSPLRKASLSSSSSLGPKERQTGAAADM
 PQPAVRHQGREPLLVKPTHSRPEGGCREVPQAASKTHGLLQAARPAQDKRPAVTSQPCPPAATHSLG
 L
 GSNLSFGPGAKRPAQAPIQACLNFPKKPRLGPFQIPESAIQGGELGAPENLQPPPAATELGSTSPQMGR
 RTPAQVPSVDRQPPHSRCLPTAQACTMSHHPAASHDGAQPLRVLFRRLENGRWSSSLLAAPSFHSP
 P
 GAFLAQSPHVSEKSEAPCVRVPPSVLYEDLQVSSSSESDSDLE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	49.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001157927</u>
Locus ID:	441315
Cytogenetics:	8p23.1
RefSeq ORF:	1392
Synonyms:	FAM90A5
Summary:	FAM90A5 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, May 2010]

Product images:



Coomassie blue staining of purified FAM90A5 protein (Cat# [TP328385]). The protein was produced from HEK293T cells transfected with FAM90A5 cDNA clone (Cat# [RC228385]) using MegaTran 2.0 (Cat# [TT210002]).