

Product datasheet for TP328522

CCDC120 (NM_001163321) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 120 (CCDC120), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228522 representing NM_001163321 Red=Cloning site Green=Tags(s)

MEVKGQLISSPTFNAPAALFGEAAPQVKSERLRGLLDRQRTLQEALSLKLQELRKVCLQEELTGQLPPE
CPLEPGERPQLVRRRPPTARAYPPPHNPQAHHSCLPAEELALEALEREVSVQQQIAAAARRLALAPDLST
EQRRRRRQVQADALRRLHELEEQLRDVRLGLPVLPLPQLPLSTGSVITTQGVCLGMRLAQLSQEDV
LHSESSLSSESASHDNEEPHGCFSLAERPSPPKAWDQLRAVSGGSPERRTPWKPPPSDLYGDLKSRNS
VASPTSPTRSLPRSASSFEGRSVPATPVLTRGAGPQLCKPEGLHSRQWSGSQDSQMGFPRADPASDRASL
FVARTRRSNSSEALLVDRAAGGGAGSPPAPLAPSASGPPVCKSSEVLYERPQTPAFSSRTAGPPDPPRA
ARPSSAAPASRGAPRLPPVCGDFLLDYSLDRGLPRSGGGTGWGELPPAAEVPGLSRRDGLLTMPLGPPPP
VYAADSNPPLLRTKDPHTRATRTKPCGLPPEAAEGPEVHPNPLLWMPPTTRIPSAGERSGHKNLALEGLR
DWYIRNSGLAAGPQRRPVLPSVGPPHPPFLHARCIEVGVQALYGAPSQAPLPHSRSTAPPVSGRYADFL
YPPELSARLSDLTLEGEQSSSDTQTPGTLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	70.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001156793
Locus ID:	90060
UniProt ID:	Q96HB5
Cytogenetics:	Xp11.23
RefSeq ORF:	1983
Synonyms:	JM11
Summary:	This gene encodes a protein that contains a coiled-coil domain. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. [provided by RefSeq, Nov 2012]

Product images:



Coomassie blue staining of purified CCDC120 protein (Cat# TP328522). The protein was produced from HEK293T cells transfected with CCDC120 cDNA clone (Cat# [RC228522]) using MegaTran 2.0 (Cat# [TT210002]).