

Product datasheet for TP328604

KLHL29 (NM_052920) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human kelch-like 29 (Drosophila) (KLHL29), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228604 representing NM_052920 Red =Cloning site Green =Tags(s)

MSRHHSRFRDYRVGWDRREWSVNGTHGTTSCSVTSGAGGGTASSLSVRPGLLPVVP SRLPTPATAP
 APCTTGSEAITSLVASSASAVTTKAPGISKGDSQSQGLATSIRWGQTPINQSTPWDTDEPPSKQMRES
 NPGTGPWVTTVAAGNQPTLIAHSYGVAQPPTFSPAVNVQAPVIGVTPSLPPHVG PQLPLMPGHYSLPQPP
 SQPLSSVVVNMPAQALYASPQPLAVSTLPGVGQVARPGPTAVGNHGMAGLLPPPPPAQPSATLPSGAPA
 TNGPPTTDSAHGLQMLRTIGVGKYEFTDPGHPREMLKELNQRRAKAFTDLKIVVEGREFEVHQNVLASC
 SLYFKDLIQRSVQDSGQGGREKLELVLSNLQADVLELLLEFVYTGSLVIDSANAKTLLEAASKFQFHTFC
 KVCVSFLEKQLTASNCLGVLAMAEAMQCSELYHMAKAFALQIFPEVAAQEEILSISKDDFIAYVSNDSLN
 TKAELVYETVIKWKKDPATRTQYAAELLAVVRLPFHPSYLLNVDNEELIKSSEACRDLVNEAKRYH
 MLPHARQEMQTPRTRPRLSAGVAEVIVLVGGRQMGVMTQSRSLVAVTCWNPQNNKWYPLASLPFYDREFFS
 VVSAGDNIYLSGGMESGVTADVWCYMSLLDNWNLVSRMTVPCRHNLSLVYDGKIYTLGGLGVAGNV DHV
 ERYDTITNQWEAVAPLPKAVHSAAATVCGGKIYVFGGVNEAGRAAGVLQSYVPQTNTWSFIESPMIDNKY
 APAVTLNGFVFILGGAYARATTIYDPEKGNKAGPNMNHRSRQFCSAVVLDGKIYATGGIVSSEGPALGNM
 EAYEPTTNTWTLLPHMPCPVFRHGCVVIKKYIQSG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	94 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_443152
Locus ID:	114818
UniProt ID:	Q96CT2
Cytogenetics:	2p24.1
RefSeq ORF:	2625
Synonyms:	KBTBD9

Product images:



Coomassie blue staining of purified KLHL29 protein (Cat# TP328604). The protein was produced from HEK293T cells transfected with KLHL29 cDNA clone (Cat# [RC228604]) using MegaTran 2.0 (Cat# [TT210002]).