

Product datasheet for **TP328604L**

KLHL29 (NM_052920) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human kelch-like 29 (Drosophila) (KLHL29), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228604 representing NM_052920 Red =Cloning site Green =Tags(s)

MSRHHSRFRDYRVGWDRREWSVNGTHGTTSCSVTSGAGGGTASSLSVRPGLPLPVVPSRLPTPATAP
APCTTGSSAITSLVASSASAVTTKAPGISKGDSQSQGLATSIRWGQTPINQSTPWDTDEPPSKQMRES
NPGTGPPWTTVAAGNQPTLIAHSYGVAQPPTFSPAVNVQAPVIGVTPSLPPHVGPQLPLMPGHYSLPQP
P
SQPLSSVVVNMPAQALYASPQPLAVSTLPGVGQVARPGPTAVGNHGMAGPLPPPPPAQPSATLPSGAP
A
TNGPPTTDSAHLQMLRTIGVGKYEFTDPGHPREMLKELNQRRAKAFTDLKIVVEGREFEVHQNVLASC
SLYFKDLIQRSVQDSGQGGREKLELVLSNLQADVLELLLEFVYTGSLVIDSANAKTLLEAASKFQHTFC
KVCVSFLEKQLTASNCLGVLAMAEAMQCSELYHMAKAFALQIFPEVAAQEEILSISKDDFIAYVSNDSLN
TKAEELVYETVIKWKIKDPATRTQYAAELLAVVRLPFIHPSYLLNVDNEELIKSSEACRDLVNEAKRYH
MLPHARQEMQTPRTRPRLSAGVAEVIVLVGGRQMVGMTQRSVLAVTCWNPQNNKWYPLASLPFYDREF
FS
VVSAGDNIYLSGGMESGVTLADVWCYMSLLDNWNLVSRMTVPRCRHNSLVYDGKIYTLGGLGVAGNV
HV
ERYDTITNQWEAVAPLPKAVHSAAATVCGGKIYVFGGVNEAGRAAGVLQSYVPQTNTWSFIESPMIDNKY
APAVTLNGFVILGGAYARATTIYDPEKGNKAGPNMNHRSRQFCSAVLDGKIYATGGIVSSEGPALGNM
EAYEPTTNTWTLLPHMPCPVFRHGCVVIKYYIQSG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	94 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_443152</u>
Locus ID:	114818
UniProt ID:	<u>Q96CT2</u>
Cytogenetics:	2p24.1
RefSeq ORF:	2625
Synonyms:	KBTBD9

Product images:



Coomassie blue staining of purified KLHL29 protein (Cat# [TP328604]). The protein was produced from HEK293T cells transfected with KLHL29 cDNA clone (Cat# [RC228604]) using MegaTran 2.0 (Cat# [TT210002]).