

Product datasheet for **TP328721M**

FXYD6 (NM_001164836) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens FXYD domain containing ion transport regulator 6 (FXYD6), transcript variant 4, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228721 protein sequence Red =Cloning site Green =Tags(s) MELVLVFLCSLLAPMVLASAAEKEKEMDPFHYDYQTLRIGGLVFAVLFVSVGILLILSRCKCSFNQKPR APGDEEAQVENLITANATEPQKAEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001158308
Locus ID:	53826
UniProt ID:	Q9H0Q3 , A0A024R3J8
RefSeq Size:	2229



[View online »](#)


Cytogenetics: 11q23.3

RefSeq ORF: 285

Summary: This gene encodes a member of the FXVD family of transmembrane proteins. This particular protein encodes phosphohippolin, which likely affects the activity of Na,K-ATPase. Multiple alternatively spliced transcript variants encoding the same protein have been described. Related pseudogenes have been identified on chromosomes 10 and X. Read-through transcripts have been observed between this locus and the downstream sodium/potassium-transporting ATPase subunit gamma (FXVD2, GeneID 486) locus.[provided by RefSeq, Feb 2011]

Protein Families: Ion Channels: Other, Transmembrane

Product images:



Coomassie blue staining of purified FXVD6 protein (Cat# [TP328721]). The protein was produced from HEK293T cells transfected with FXVD6 cDNA clone (Cat# [RC228721]) using MegaTran 2.0 (Cat# [TT210002]).