

## Product datasheet for **TP328740L**

### **FHIT (NM\_001166243) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fragile histidine triad gene (FHIT), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228740 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MSFRFGQHLLIKPSVFLKTELSFALVNRKPVVPGHVLVCLRPVERFHDLRPDEVADLFQTTQRVGTVE KHFHGTSLTFSMQDGPEAGQTVKHVHVHVLPRKAGDFHRNDSIYEELQKHKEDFPASWRSEEEEMAAEAA ALRVYFQ  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	16.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001159715</a>
Locus ID:	2272
UniProt ID:	<a href="#">P49789</a> , <a href="#">A0A024R366</a>
RefSeq Size:	1122



[View online »](#)

Cytogenetics: 3p14.2


RefSeq ORF: 441

Synonyms: AP3Aase; FRA3B

**Summary:** The protein encoded by this gene is a P1-P3-bis(5'-adenosyl) triphosphate hydrolase involved in purine metabolism. This gene encompasses the common fragile site FRA3B on chromosome 3, where carcinogen-induced damage can lead to translocations and aberrant transcripts. In fact, aberrant transcripts from this gene have been found in about half of all esophageal, stomach, and colon carcinomas. The encoded protein is also a tumor suppressor, as loss of its activity results in replication stress and DNA damage. [provided by RefSeq, Aug 2017]

**Protein Pathways:** Non-small cell lung cancer, Purine metabolism, Small cell lung cancer

### Product images:



Coomassie blue staining of purified FHIT protein (Cat# [TP328740]). The protein was produced from HEK293T cells transfected with FHIT cDNA clone (Cat# [RC228740]) using MegaTran 2.0 (Cat# [TT210002]).