

Product datasheet for **TP328766**

C21orf62 (NM_001162495) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens chromosome 21 open reading frame 62 (C21orf62), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228766 representing NM_001162495 Red =Cloning site Green =Tags(s)
	 MAPPSRHCLLLISTLGVFALNCFTKGQKNSTLIFTRENTIRNCSCSADIRDCDYSLANLMCNCKTVLPLA VERTSYNGHLTIWFTDTSALGHLLNFTLVQDLKLSLCSTNTLPTEYLAICGLKRLRINMEAKHPFPEQSL LIHSGGSDSREKPMWLHKGWQPCMYISFLDMALFNDRSALKSYSIENVTSIANNFPDFSYFRTPMPSN KSYWTFIY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001155967</u>
Locus ID:	56245



[View online »](#)

UniProt ID: [Q9NYP8](#)
Cytogenetics: 21q22.11
RefSeq ORF: 657
Synonyms: B37; C21orf120; PRED81

Product images:



Coomassie blue staining of purified C21orf62 protein (Cat# TP328766). The protein was produced from HEK293T cells transfected with C21orf62 cDNA clone (Cat# [RC228766]) using MegaTran 2.0 (Cat# [TT210002]).