

Product datasheet for TP328855

HSD3B2 (NM_001166120) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 2 (HSD3B2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC228855 protein sequence Red =Cloning site Green =Tags(s) MGWSCLVTGAGGLLGQRIVRLVVEEKELKEIRALDKAFRPELREEFSKLQNRTKLTVLEGDILDEPFLKR ACQDVSIVIHTACIIDVFGVTHRESIMNVNVKGTQLLLEACVQASVPVFIYTSSIEVAGPNSYKEIIQNG HEEPLENTWPTPYPSKKLAEKAVLAANGWNLKNGDTLYTCALRPTYIYGEGGPFLSASINEALNNNGI LSSVGKFSTVNPVYVGNVAWAHILALRALRDPKKAPSVRGQFYIISDDTPHQSYDNLNLYLSKEFGLRLD SRWSLPLTLMYWIGFLLEVVSFLLSPIYSYQPPFNRHTVTLNSVFTFSYKKAQRDLAYKPLYSWEEAKQ KTVWVGSLVDRHKETLKSQTQ TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	41.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001159592</u>


[View online »](#)

Locus ID:	3284
UniProt ID:	<u>P26439</u>
RefSeq Size:	1807
Cytogenetics:	1p12
RefSeq ORF:	1116
Synonyms:	HSD3B; HSDB; SDR11E2
Summary:	The protein encoded by this gene is a bifunctional enzyme that catalyzes the oxidative conversion of delta(5)-ene-3-beta-hydroxy steroid, and the oxidative conversion of ketosteroids. It plays a crucial role in the biosynthesis of all classes of hormonal steroids. This gene is predominantly expressed in the adrenals and the gonads. Mutations in this gene are associated with 3-beta-hydroxysteroid dehydrogenase, type II, deficiency. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Oct 2009]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Androgen and estrogen metabolism, C21-Steroid hormone metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified HSD3B2 protein (Cat# TP328855). The protein was produced from HEK293T cells transfected with HSD3B2 cDNA clone (Cat# [RC228855]) using MegaTran 2.0 (Cat# [TT210002]).