

## Product datasheet for TP329090L

### TIMM8B (NM\_012459) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human translocase of inner mitochondrial membrane 8 homolog B (yeast) (TIMM8B), nuclear gene encoding mitochondrial protein, transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229090 representing NM_012459 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MRKHSCRKVASLRRTMAELGEADEAELQRLVAAEQQKAQFTAQVHHFMELCWDKCVEKPGNRLDSRTE NC LSSCVDRFIDTTLAITSRFAQIVQKGGQ  <b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	11 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_036591</a></u>
Locus ID:	26521
UniProt ID:	<u><a href="#">Q9Y5J9</a></u>


[View online »](#)


**Cytogenetics:** 11q23.1

**RefSeq ORF:** 294

**Synonyms:** DDP2; TIM8B

**Summary:** This gene encodes a member of a well-conserved family of proteins with similarity to yeast Tim mitochondrial import proteins. This gene is encoded by a nuclear gene and is transported into the intermembrane space of the mitochondrion. When formed into complexes, these proteins guide membrane-spanning proteins across the mitochondrial intermembrane space before they are added into the mitochondrial inner membrane. This gene is adjacent to succinate dehydrogenase, subunit D (SDHD), in which mutations have been found in affected members of families with hereditary paraganglioma.[provided by RefSeq, Aug 2009]

### Product images:



Coomassie blue staining of purified TIMM8B protein (Cat# [TP329090]). The protein was produced from HEK293T cells transfected with TIMM8B cDNA clone (Cat# [RC229090]) using MegaTran 2.0 (Cat# [TT210002]).