

## Product datasheet for TP329117

### Interferon alpha (IFNA13) (NM\_006900) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens interferon, alpha 13 (IFNA13), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229117 representing NM_006900 Red=Cloning site Green=Tags(s)
	MASPFALLMALVVLVLSCKSSCSLGCGLPETHSLDNRRTLMLLAQMSRISPSSCLMDRHDGFPQEEFDGNQ FQKAPAI SVLHELIIQQIFNLFTTKDSSAAWDEDLLDKFCTELYQQLNDLEACVMQEERVGETPLMNADSI LAVKKYFRITLYLTEKKYSPCAWEVVRAEIMRSLSLSTNLQERLRKE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_008831</a>
Locus ID:	3447
UniProt ID:	<a href="#">P01562</a> , <a href="#">A0A087WWS6</a>
Cytogenetics:	9p21.3



[View online »](#)


RefSeq ORF: 567

**Summary:** Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. [UniProtKB/Swiss-Prot Function]

**Protein Families:** Druggable Genome, Secreted Protein

**Protein Pathways:** Antigen processing and presentation, Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Cytosolic DNA-sensing pathway, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of autophagy, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway

### Product images:



Coomassie blue staining of purified IFNA13 protein (Cat# TP329117). The protein was produced from HEK293T cells transfected with IFNA13 cDNA clone (Cat# [RC229117]) using MegaTran 2.0 (Cat# [TT210002]).