

Product datasheet for TP329124

LMO2 (NM_005574) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human LIM domain only 2 (rhombotin-like 1) (LMO2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229124 representing NM_005574 Red =Cloning site Green =Tags(s)
	<p>MEGSAVTVLERGGASSPAERRSKRRRRSGDGGGGGGARAPEGVRAPAAGQPRATKGAPPPPGTPPPSPM SSAIERKSLDPSEEPVDEVLQIPPSLLTCGGCQQNIGDRYFLK AIDQYWHEDCLSCDLGCGRLGEVGRRL YYKLGRKLCRRDYLR LFGQDGLCASCDKRIRAYEMTMRVKDKVYHLECFKCAACQKHFCVGDYLLINS IVCEQDIYEWTKINGMI</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	24.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_005565</u>
Locus ID:	4005



[View online »](#)

UniProt ID:	P25791
Cytogenetics:	11p13
RefSeq ORF:	681
Synonyms:	LMO-2; RBTN2; RBTNL1; RHOM2; TTG2
Summary:	LMO2 encodes a cysteine-rich, two LIM-domain protein that is required for yolk sac erythropoiesis. The LMO2 protein has a central and crucial role in hematopoietic development and is highly conserved. The LMO2 transcription start site is located approximately 25 kb downstream from the 11p13 T-cell translocation cluster (11p13 ttc), where a number T-cell acute lymphoblastic leukemia-specific translocations occur. Alternative splicing results in multiple transcript variants encoding different isoforms.[provided by RefSeq, Nov 2008]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified LMO2 protein (Cat# TP329124). The protein was produced from HEK293T cells transfected with LMO2 cDNA clone (Cat# [RC229124]) using MegaTran 2.0 (Cat# [TT210002]).