

## Product datasheet for TP329177

### PRMT6 (NM\_018137) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein arginine methyltransferase 6 (PRMT6), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229177 representing NM_018137 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSQPKKRKLESGGGEGGEGTEEDGAEREAALERPRRTKRERDQLYYECYSDVSVHEEMIADRVRTDAY RLGILRNWAALRGKTVLDVGAGTGILSIFCAQAGARRVYAVEASAIWQQAREVVRFNGLEDRVHVLPGPV ETVELPEQVDAIVSEWMGYLLHESMLSSVLHARTKWLKEGGLLLPASAELFIAPISDQMLEWRLGFWSQ VKQHYGVDMSCLEGFATRCLMGHSEIVVQGLSGEDVLRPQRFAQLELSRAGLEQELEAGVGGRFCSCY GSAPMHGFAIWFQVTFPGGESEKPLVLSTSPFHPATHWKQALLYLNEPVQVEQDQDVSGETLLPSRDNP RRLRVLLRYKVGQDQEEKTKDFAMED
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	41.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_060607</a>
Locus ID:	55170



[View online »](#)

UniProt ID: [Q96LA8](#)

Cytogenetics: 1p13.3


RefSeq ORF: 1125

Synonyms: HRMT1L6

**Summary:** The protein encoded by this gene belongs to the arginine N-methyltransferase family, which catalyze the sequential transfer of methyl group from S-adenosyl-L-methionine to the side chain nitrogens of arginine residues within proteins, to form methylated arginine derivatives and S-adenosyl-L-homocysteine. This protein can catalyze both, the formation of omega-N monomethylarginine and asymmetrical dimethylarginine, with a strong preference for the latter. It specifically mediates the asymmetric dimethylation of Arg2 of histone H3, and the methylated form represents a specific tag for epigenetic transcriptional repression. This protein also forms a complex with, and methylates DNA polymerase beta, resulting in stimulation of polymerase activity by enhancing DNA binding and processivity. [provided by RefSeq, Sep 2011]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified PRMT6 protein (Cat# TP329177). The protein was produced from HEK293T cells transfected with PRMT6 cDNA clone (Cat# [RC229177]) using MegaTran 2.0 (Cat# [TT210002]).