

Product datasheet for **TP329177L**

PRMT6 (NM_018137) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human protein arginine methyltransferase 6 (PRMT6), 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC229177 representing NM_018137
Red=Cloning site **Green**=Tags(s)

MSQPKKRKLESGGGEGGEGTEEDGAEREAALERPRRTKRERDQLYYECYSDVSVHEEMIADRVRTDAY
RLGILRNWAALRGKTVLDVGAGTGILSIFCAQAGARRVYAVEASAIWQQAREVVRFNGLDRVHVLPGPV
ETVELPEQVDAIVSEWMGYGLLHESMLSSVLHARTKWLKEGGLLPASAELFIAPISDQMLEWRLGFWSQ
VKQHYGVDMSCLEGFATRCLMGHSEIVVQGLSGEDVLRPQRFAQLELSRAGLEQELEAGVGGFRFCSCY
GSAPMHGFAIWFQVTFPGGESEKPLVLSTSPFHPATHWKQALLYLNEPVQVEQDQTDVSGEITLLPSRDNP
RRLRVLLRYKVGQDQEETKDFAMED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 41.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_060607](#)

Locus ID: 55170



[View online »](#)

UniProt ID: [Q96LA8](#)

Cytogenetics: 1p13.3


RefSeq ORF: 1125

Synonyms: HRMT1L6

Summary: The protein encoded by this gene belongs to the arginine N-methyltransferase family, which catalyze the sequential transfer of methyl group from S-adenosyl-L-methionine to the side chain nitrogens of arginine residues within proteins, to form methylated arginine derivatives and S-adenosyl-L-homocysteine. This protein can catalyze both, the formation of omega-N monomethylarginine and asymmetrical dimethylarginine, with a strong preference for the latter. It specifically mediates the asymmetric dimethylation of Arg2 of histone H3, and the methylated form represents a specific tag for epigenetic transcriptional repression. This protein also forms a complex with, and methylates DNA polymerase beta, resulting in stimulation of polymerase activity by enhancing DNA binding and processivity. [provided by RefSeq, Sep 2011]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified PRMT6 protein (Cat# [TP329177]). The protein was produced from HEK293T cells transfected with PRMT6 cDNA clone (Cat# [RC229177]) using MegaTran 2.0 (Cat# [TT210002]).