

Product datasheet for TP329260

UAP1L1 (NM_207309) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human UDP-N-acetylglucosamine pyrophosphorylase 1-like 1 (UAP1L1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229260 representing NM_207309 Red =Cloning site Green =Tags(s)

MASEQDVRARLQRAGQEHLRFWAELAPEPRAALLAELALLEPEALREHCRRAAEACARPHGPPDLAAR
LRPLPPERVGRASRSDPETRRRWEEEGFRQISLNKVAVLLLAGGQGTRLGVTPKGMVYRVLPSRKTLYQ
LQAERIRRVEQLAGERHGTCTVPWYVMTSEFTLGPTAEFFREHNFFHLDPANVVMFEQRLLPVAVTFDGK
VILERKDKVAMAPDGNGLYCALEDHKILEDMERRGVEFVHVYCVDNILVRLADPVFIGFVQLQGADCGA
KVVEKAYPEEPVGVWCQVDGVPQVVEYSEISPETAQLRASDGSLLYNAGNICNHFFTRGFLKAVTREFEP
LLKPHVAVKKVPYVDEEGLVKPLKPNGIKMEKFVDFVFRFAKNFAALEVLRREEFSPLKNAEPADRSDP
RTARQALLTQHYRWALRAGARFLDAHGAWLPELPSLPPNGDPPAICEISPLVSYSGEGLEVYLQGREFQS
PLILDEDQAREPQLQES


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	56.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_997192
Locus ID:	91373
UniProt ID:	Q3KQV9
Cytogenetics:	9q34.3
RefSeq ORF:	1521

Product images:

Coomassie blue staining of purified UAP1L1 protein (Cat# TP329260). The protein was produced from HEK293T cells transfected with UAP1L1 cDNA clone (Cat# [RC229260]) using MegaTran 2.0 (Cat# [TT210002]).