

Product datasheet for **TP329534L**

C6orf223 (NM_001171992) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 223 (C6orf223), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229534 representing NM_001171992 Red =Cloning site Green =Tags(s) MMPLAEAGALAQGGGPSATEWACILRRLRLKEGVDGGFLFSPEDPQAQAAHLTDGQSIQKKWQDFSCPQSR KTECFWQKEQHKQRSGNPRCLQFEGGERTPPPPQT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	12.1
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001165463</u>
Locus ID:	221416
Cytogenetics:	6p21.1




[View online »](#)

RefSeq ORF: 315

Synonyms: MGC45491

Product images:



Coomassie blue staining of purified C6orf223 protein (Cat# [TP329534]). The protein was produced from HEK293T cells transfected with C6orf223 cDNA clone (Cat# [RC229534]) using MegaTran 2.0 (Cat# [TT210002]).