


Product datasheet for **TP329555M**

FAM229A (NM_001167676) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical LOC100128071 (LOC100128071), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229555 representing NM_001167676 Red =Cloning site Green =Tags(s) MLPSSTPGPGHATETCPAPPGPERSPAARAPAAASSLGPVSTAGRAPRGLDMSAQEPPQGRRFPPIEAGDS RGLAAPESQDSPEAVATEHNVPVRPLRRCPGCHCLTLLHVPIDVYLAMGGSPRARAT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.4
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001161148
Locus ID:	100128071
UniProt ID:	H3BQW9
Cytogenetics:	1p35.1
RefSeq ORF:	381


[View online »](#)

Product images:

Coomassie blue staining of purified FAM229A protein (Cat# [TP329555]). The protein was produced from HEK293T cells transfected with FAM229A cDNA clone (Cat# [RC229555]) using MegaTran 2.0 (Cat# [TT210002]).