

Product datasheet for **TP329599M**

Aquaporin 1 (AQP1) (NM_001185062) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens aquaporin 1 (Colton blood group) (AQP1), transcript variant 4, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229599 representing NM_001185062 Red =Cloning site Green =Tags(s)
	MQSGMGWNVLDLFWLADGVNSGQGLGIEIIGTLQLVLCVLATDRRRRDLGGSAPLAIGLSVALGHLLAID YTGCGINPARSFGSAVITHNFSNHWFVWGPFIGGALAVLIYDFILAPRSSDLTDRVKVWTSQVEEYDL DADDINSRVEMKPK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.1
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001171991</u>
Locus ID:	358
UniProt ID:	<u>P29972</u>



[View online »](#)

Cytogenetics: 7p14.3


RefSeq ORF: 462

Synonyms: AQP-CHIP; CHIP28; CO

Summary: This gene encodes a small integral membrane protein with six bilayer spanning domains that functions as a water channel protein. This protein permits passive transport of water along an osmotic gradient. This gene is a possible candidate for disorders involving imbalance in ocular fluid movement. [provided by RefSeq, Aug 2016]

Protein Families: Druggable Genome, Ion Channels: Other, Transmembrane

Product images:



Coomassie blue staining of purified AQP1 protein (Cat# [TP329599]). The protein was produced from HEK293T cells transfected with AQP1 cDNA clone (Cat# [RC229599]) using MegaTran 2.0 (Cat# [TT210002]).