

Product datasheet for TP329617L

PQBP1 (NM_001167992) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens polyglutamine binding protein 1 (PQBP1), transcript variant 10, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229617 representing NM_001167992 Red =Cloning site Green =Tags(s) MPLPVALQTRLAKRGILKHLEPEPEEEIIAEDYDDDPVDYEATRLEGLPPSWYKVFDPSCGLPYYWNADR EEGKERRHHRREELAPYPKSKKAVSRKDEELDPMDPSSYSADAPRGTWSTGLPKRNEAKTGADTTAAGPLF QQRYPSPGAVLRANAASRTKQQD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.2
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001161464</u>
Locus ID:	10084
UniProt ID:	<u>O60828</u>



[View online »](#)

Cytogenetics: Xp11.23

RefSeq ORF: 495


Synonyms: MRX2; MRX55; MRXS3; MRXS8; NPW38; RENS1; SHS

Summary: This gene encodes a nuclear polyglutamine-binding protein that is involved with transcription activation. The encoded protein contains a WW domain. Mutations in this gene have been found in patients with Renpenning syndrome 1 and other syndromes with X-linked cognitive disability. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.[provided by RefSeq, Nov 2009]

Protein Families: Transcription Factors

Protein Pathways: Spliceosome

Product images:



Coomassie blue staining of purified PQBP1 protein (Cat# [TP329617]). The protein was produced from HEK293T cells transfected with PQBP1 cDNA clone (Cat# [RC229617]) using MegaTran 2.0 (Cat# [TT210002]).