

Product datasheet for **TP329654**

TSPAN3 (NM_001168412) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens tetraspanin 3 (TSPAN3), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229654 representing NM_001168412 Red =Cloning site Green =Tags(s)
	 MGQCGITSSKTVLVFLNLIWFVVIILLVVFTEVWWWLGYVYRAKVENEVDRSIQKVYKTYNGTNPDAASRAIDYVQRQLHCCGIHNYSDWENTDWFKETKNQSVPLSCCRETASNCNGSLAHPSDLYAEGCEALVWKKLQEIMMHVIWAALAFAAIQLLGMCLCACIVLCRRSRDPAYELLITGGTYA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.6
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001161884</u>
Locus ID:	10099
UniProt ID:	<u>O60637</u>



[View online »](#)

Cytogenetics: 15q24.3


RefSeq ORF: 567

Synonyms: TM4-A; TM4SF8; TSPAN-3

Summary: The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The use of alternate polyadenylation sites has been found for this gene. Multiple alternative transcripts encoding different isoforms have been described. [provided by RefSeq, Dec 2009]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified TSPAN3 protein (Cat# TP329654). The protein was produced from HEK293T cells transfected with TSPAN3 cDNA clone (Cat# [RC229654]) using MegaTran 2.0 (Cat# [TT210002]).