

## Product datasheet for TP329656L

### CD99L2 (NM\_001184808) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens CD99 molecule-like 2 (CD99L2), transcript variant 4, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229656 representing NM_001184808 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MVAWRS AFLVCLAFSLATLVQRGSGDFDDFNLEDAVKETSSVKQPWDHTTTTTTNRPGTTRAPAKPPGSG LDLADALDDQDDGRRKPGIGGRGDGRYGSNDDPGSGMVAEPGTIAGVASALAMALIGAVSSYISYQQKF CFSIQQLNADYVKGENLEAVVCEEPQVKYSTLHTQSAEPPPPPEPARI  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	20.5
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001171737</a>
Locus ID:	83692
UniProt ID:	<a href="#">Q8TCZ2</a>



[View online »](#)

Cytogenetics: Xq28


RefSeq ORF: 567

Synonyms: CD99B; MIC2L1

**Summary:** This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]

**Protein Families:** Transmembrane

### Product images:



Coomassie blue staining of purified CD99L2 protein (Cat# [TP329656]). The protein was produced from HEK293T cells transfected with CD99L2 cDNA clone (Cat# [RC229656]) using MegaTran 2.0 (Cat# [TT210002]).