

## Product datasheet for TP329759

### C3orf80 (NM\_001168214) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC401097 (LOC401097), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229759 representing NM_001168214 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MWGPGVTAEGLSVAPAPPPLLPLLLLLLALALVAPSRGGGGCAELACGERERCCDATNATAVRCKLPLHA FLDNVGFVVRKLSGLLILLVLFAIGYFLQRIICPSPRRYPRGQARPGQRPGGAGPLGGAGPPDDDDD SPALLRDEAAAGSQDSLDSGGGGRGRGGGGRSDPSCASEHEMRVSPVFLQLPSYEEVKYLPTYEESMR LQQLSPGEVLPVSVLGRPRGGVAAEPDGGEGRYPLI  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	26.1
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001161686</a>
Locus ID:	401097
UniProt ID:	<a href="#">F5H4A9</a>




[View online »](#)

Cytogenetics: 3q25.33

RefSeq ORF: 741

Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified C3orf80 protein (Cat# TP329759). The protein was produced from HEK293T cells transfected with C3orf80 cDNA clone (Cat# [RC229759]) using MegaTran 2.0 (Cat# [TT210002]).