

Product datasheet for **TP329840**

PGAM5 (NM_001170543) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphoglycerate mutase family member 5 (PGAM5), nuclear gene encoding mitochondrial protein, transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229840 representing NM_001170543 Red =Cloning site Green =Tags(s)
	<p>MAFRQALQLAACGLAGGSAAVLFSAVAVGKPRAGGDAEPRPAEPPAWAGGARPGPGVWDPNWRREPLSL INVRKRNVESGEEELASKLDHYKAKATRHIFLIRHSQYHVDGSLEKDRTLPLGREQAELTGLRLASLGL KFNKIVHSSMTRAIETTDIISRHLPGVCKVSTDLLREGAPIEPDPPVSHWKPEAVQYEDGARIEAAFRN YIHRADARQEEDSYEIFICHANVIRYIVCRALQFPPEGWLRLSLNNGSITHLVIRPNGRVALRTLGDGTGF MPPDKITRS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	32.5
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001164014</u>
Locus ID:	192111



[View online »](#)

UniProt ID: [Q96HS1](#)

Cytogenetics: 12q24.33


RefSeq ORF: 867

Synonyms: BXLBV68

Summary: Displays phosphatase activity for serine/threonine residues, and, dephosphorylates and activates MAP3K5 kinase. Has apparently no phosphoglycerate mutase activity. May be regulator of mitochondrial dynamics. Substrate for a KEAP1-dependent ubiquitin ligase complex. Contributes to the repression of NFE2L2-dependent gene expression. Acts as a central mediator for programmed necrosis induced by TNF, by reactive oxygen species and by calcium ionophore. [UniProtKB/Swiss-Prot Function]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified PGAM5 protein (Cat# TP329840). The protein was produced from HEK293T cells transfected with PGAM5 cDNA clone (Cat# [RC229840]) using MegaTran 2.0 (Cat# [TT210002]).