

Product datasheet for **TP329854L**

MBNL3 (NM_001170701) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens muscleblind-like 3 (Drosophila) (MBNL3), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229854 representing NM_001170701 Red =Cloning site Green =Tags(s) MERASKNLKGRCTRENCKYLHPPPHLKTQLEINGRNLIQQKTAAMFAQQMQLMLQNAQMSSLGSPMT PSIPANPPMAFNPIPHPGMGLVPAELVPNTVLPVLPGNPPLAMPGAVGPKLMRSDKLEVCREFQRGNCTR GENDCRYAHPTDASMIEASDNTVTICMDYIKGRCSREKCKYFHPPAHLQARLKAHHQMNHSAASAMALQ PGTLQLIPKRSALPKNGATPVFNPTVFHCQQALTNLQLPQPAFIPAGPILCMAPASNIVPMMHGATPTT VSAATTPATSVFPAAPTGNQLKF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.4
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001164172</u>
Locus ID:	55796



[View online »](#)

UniProt ID: [Q9NUK0](#)


Cytogenetics: Xq26.2

RefSeq ORF: 912

Synonyms: CHCR; MBLX; MBLX39; MBXL

Summary: This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described. [provided by RefSeq, Dec 2009]

Product images:



Coomassie blue staining of purified MBNL3 protein (Cat# [TP329854]). The protein was produced from HEK293T cells transfected with MBNL3 cDNA clone (Cat# [RC229854]) using MegaTran 2.0 (Cat# [TT210002]).