

Product datasheet for TP329902L

CYB5R3 (NM_001171660) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens cytochrome b5 reductase 3 (CYB5R3), transcript variant 5, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC229902 representing NM_001171660 Red =Cloning site Green =Tags(s) MNRSLLVGCMQSKDIWGREESICERLKQDGLDVERAESWELGHMVLPVWFLYLLMKLFQRSTPAITL SPDIKYPLRLIDREIISHDTRRRFALPSPQHILGLPVGQHIYLSARIDGNLVVRPYTPISSDDDKGFVD LVIKVYFKDTHPKFPAGGKMSQYLESMQIGDTIEFRGPSGLLVYQGKGKFAIRPDKKSNIIRTVKSVGM IAGGTGITPMLQVIRAIMKDPDDHTVCHLLFANQTEKDILLRPELEELRNKHSARFKLWYTLDRAPAWD YGQGFVNEEMIRDHLPPEEEPLVLMCGPPPMIQYACLPNLDHVGHPTERCFVF TR TRPLEQKLISEEDLA ANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	38.7
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001165131</u>
Locus ID:	1727


[View online »](#)

UniProt ID: P00387

Cytogenetics: 22q13.2

RefSeq ORF: 1002


Synonyms: B5R; DIA1

Summary: This gene encodes cytochrome b5 reductase, which includes a membrane-bound form in somatic cells (anchored in the endoplasmic reticulum, mitochondrial and other membranes) and a soluble form in erythrocytes. The membrane-bound form exists mainly on the cytoplasmic side of the endoplasmic reticulum and functions in desaturation and elongation of fatty acids, in cholesterol biosynthesis, and in drug metabolism. The erythrocyte form is located in a soluble fraction of circulating erythrocytes and is involved in methemoglobin reduction. The membrane-bound form has both membrane-binding and catalytic domains, while the soluble form has only the catalytic domain. Alternate splicing results in multiple transcript variants. Mutations in this gene cause methemoglobinemias. [provided by RefSeq, Jan 2010]

Protein Families: Druggable Genome

Protein Pathways: Amino sugar and nucleotide sugar metabolism

Product images:



Coomassie blue staining of purified CYB5R3 protein (Cat# [TP329902]). The protein was produced from HEK293T cells transfected with CYB5R3 cDNA clone (Cat# [RC229902]) using MegaTran 2.0 (Cat# [TT210002]).