

## Product datasheet for **TP330029L**

### ATXN10 (NM\_001167621) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ataxin 10 (ATXN10), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC230029 representing NM_001167621 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAAPRPPPARLSGVMVPAPIQDLEALRALTALFKEQRNRNLDITIGVAVDLILLFREL RVEQESLLTAFRC GLQFLGNIASRNEDSQSIVVWHAFPELFLSCLNHPDKKIVAYSSMILFTSLNHERMKELEENLNIAIDVI DAYQKHPESEWPFLIITDLFLKSPQLVQAMFPKLNQERVTLDDLMIAKITSDEPLTKDDIPVFLRHAEL IASTFVDQCKTVLKLASEPPDDEEALATIRLLDVLCEMTVNTTELLGYLQVFPGLLERVIDLLRVIHVAG KETTNI FSNCGCVRAEGDISNVANGFKSHLIRLIGNLCYKNKDNQDKVNELDGIPLILDNCNISDSNPFL TQWVIYAIRNLTEDNSQNQDLIAKMEEQGLADASLLKKVGFVEKKGEKLIKSTRDTPKP</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	46.7
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001161093</a>
Locus ID:	25814



[View online »](#)

UniProt ID: [Q9UBB4, A6NLC4](#)


Cytogenetics: 22q13.31

RefSeq ORF: 1233

Synonyms: E46L; HUMEEP; SCA10

**Summary:** This gene encodes a protein that may function in neuron survival, neuron differentiation, and neuritogenesis. These roles may be carried out via activation of the mitogen-activated protein kinase cascade. Expansion of an ATTCT repeat from 9-32 copies to 800-4500 copies in an intronic region of this locus has been associated with spinocerebellar ataxia, type 10. Alternatively spliced transcript variants have been described.[provided by RefSeq, Jul 2016]

### Product images:



Coomassie blue staining of purified ATXN10 protein (Cat# [TP330029]). The protein was produced from HEK293T cells transfected with ATXN10 cDNA clone (Cat# [RC230029]) using MegaTran 2.0 (Cat# [TT210002]).