

Product datasheet for TP330517L

C1orf127 (NM_001170754) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 127 (C1orf127), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC230517 representing NM_001170754 Red =Cloning site Green =Tags(s)

MWGGSPALAWAVWLACVQPTVFPWLSFRSDTDKPSSAAEVLTEASSCWTDTVECFSDYMTLWIPRSHVE
G
LRRWLARTLHLPGTWRSPDHL DSSLAKCGYFLHLASDGDFLFRVQYSACFVQKEKANYRLEIRIFQKGV
MLERSDRYIMKCPMLRSRLGQESVHCGPMFIQVSRPLPLWRDNRQTPWLLSLRGELVASLEDASLMGLY
V
DMNATTVTVQSPRQGLLQRWEVLNTSAELLPLWLVS GHHAYSLEAACPPVVSFQPESEVLVHIPKQRLGLV
KRGSYIEETLSRFLRVHQSNIFMV TENKDFVVSIPAAGVLQVQRCQEVGGTPTGQAFYRVDLSLEFAE
MAAPVLWTVESFFQCVGSGTESPASTAALRTTSPSPGPETPPAGVPPAASSQVWAAGPAAQEWLSRDL
LHRPSDALAKKGLGPFLQTAKPARRGQTSASILPRV VQAQRGPQPPPGEAGIPGHPTPPATLPSEPVEGV
QASPWRPRPVLPTHALTLPVSSDASSPSPAPRPERPE SLLVSGPSVTLTEGLGTVRPEQDPAKSPGSP
LLLRGLSSGDVAPEPIMGEPGQASEEFQPLARPWRATLAAEELVSHRSPGEPQETCSGTEVERPRQTGP
GLPREGARGHMDLSSSEPSQDIEGGLSILPARDATFSTPSVRQPDPSAWLSSGPELTGMPRVRLAAPLA
VLPMEPLPPEPVRPAALLTPEASSVGGPDQARYLESAPGWVPGQE EWGVAHTSSPPSTQTLSLWAPTGV
LPSLVELEYPFQAGRGASLQQELTEPTLALS AESHRPPELQDSVEGLSERPSR

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	89.5
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001164225](#)


Locus ID: 148345

UniProt ID: [G8JLG8](#)

Cytogenetics: 1p36.22

RefSeq ORF: 2469

Product images:



Coomassie blue staining of purified C1orf127 protein (Cat# [TP330517]). The protein was produced from HEK293T cells transfected with C1orf127 cDNA clone (Cat# [RC230517]) using MegaTran 2.0 (Cat# [TT210002]).