

Product datasheet for TP330517M

C1orf127 (NM_001170754) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 127 (C1orf127), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC230517 representing NM_001170754 Red=Cloning site Green=Tags(s)

MWGSPALAWAVWLACVQPTVFPWLSFRSDTDKPSAAEVLTEASSCWTDTVECFSDYMTLWIPRSHVEG
LRRWLARTLHLPGTWRSPDHL DSSLAKCGYFLHLASDGDFLFRVQYSACFVQKEKANYRLEIRIFQKGV
GLERSDRYIMKCPMLRSRLGQESVHCGPMFIQVSRPLPLWRDNRQTPWLLSLRGELVASLEDASLMGLYV
DMNATTVTVQSPRQGLLQRWEVLNTSAELLPLWLVS GHHAYSLEAACPPVSFQPESEVLVHIPKQRLGLV
KRGSYIEETLSRFLRVHQSNI FMTENKDFVVSIPAAGVLQVQRCQEVGGTPTGTA FYRVDLSLEFAE
MAAPVLWTVESFFQCVGSGTESPASTAALRTTPSPSPGPETPPAGVPPAASSQVWAAGPAAQEWLSRDL
LHRPSDALAKKGLGPFLQTAKPARRGQTSASILPRVQAQRGPQPPG EAGIPGHPTPPATLPSEVVEGV
QASPWRPRPVLPTHAL TLPVSSDASSPSPAPRPERPE LLSVSGPSVTLTEGLGTVRPEQDPAKSPGSP
LLLRLSSGDVA APEPIMGEPGQASEEFQPLARPWRATLAAEELVSHRSPGEPQETCSGTEVERPRQTGP
GLPREGARGHMDLSSSEPSQDIEGPGLSILPARDATFSTPSVRQPDPSAWLSSGPELTGM PRVRLAAPLA
VLPMEPLPPEPVRPAALLTPEASSVGGPDQARYLESAPGWPVQEEWGWVAHTSSPSTQTL SLWAPTGV
LPSLVELEYPFQAGRGASLQQLTEPTLALSAESHRPPELQDSVEGLSERPSR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	89.5
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001164225
Locus ID:	148345
UniProt ID:	G8JLG8 , B7ZLG7
Cytogenetics:	1p36.22
RefSeq ORF:	2469

Product images:



Coomassie blue staining of purified C1orf127 protein (Cat# [TP330517]). The protein was produced from HEK293T cells transfected with C1orf127 cDNA clone (Cat# [RC230517]) using MegaTran 2.0 (Cat# [TT210002]).