

## Product datasheet for TP330681L

### FAM21C (WASHC2C) (NM\_001169107) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens family with sequence similarity 21, member C (FAM21C), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC230681 representing NM_001169107 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MMNRTTPDQELVPASEPVWERPWSVEEIRRSSQSWSLAADAGLLQFLQEFSSQQTISRTHEIKKQVDGLIR  
ETKATDCRLHNVFNDFLMLSNTQFIENRVYDEEVEEVLKAEAEKTEQEKTREQKEVDLIPKVQEAENVY  
LQVLDSAFEQLDIKAGNSDSEEDDANGRVELILEPKDLYDRPLPYLIGSKLFMEQEDVGLGELSSEEGS  
VGSDRGSIVDTEEEKEEEDDFAHHSNEQNQHTTQMSDEEEDDDGCDLFADSEKEEEDIEDIEENTR  
PKRSRPTFADELAARIKGDAMGRVDEEPTNEEDNLFAPPKLTDEDFSPFGSGGGLFSGGKGLFDDEDEE  
SDLFTEASQDRQAGASVKEESSSSKPGKKIPAGAVSVFLGDTDFVGAASVPSLKEPQKPEQPTPRKSPYG  
PPPTGLFDDDDGDDDDFFSAPHSKPSKTRKVQSTADIFGDEEGDLFKEKAVASPEATVSQTDENKARAE  
KKVTLSSYKLNKPSSETKTQKGLFSDEEEDSEDLFSSQSASNLKGASLLPGKLPVSFLDDEDEEDNLF  
GTAACKQTLQLAQREEKAKASELSKKKASALLFSSDEEDQWNIPASQTHLASDSRSKGEPRDSGTLQSQ  
EAKAVKKTSLFEEDKEDDLFAIAKDSQKKTQRVSLLFEDDVSGLSFGSPPTSVPATKEGLLTRSAQE  
TVKHSDLFSSSSPWDKGTKPRTKTVLSLFDEEEDKMEDQNIIQAPQKEVKGKCDPDAHPKSTGVFQDEEL  
LFSHKLQKDNPDVDFLAFGTTKTKLLEPSVGLSFGDDEDDDLFSSAKSQPLVQEKKRVVKKDHSVNSFKN  
QKHPESIQGSKEKGIWKPETPQANLAINPAALLPTAASQISEVKPVLPELAFPSSEHRRSHGLESVPVLP  
GSGEAGVSFDLPAQADTLHSANKSRVKMRGKRRPQTRAARRLAAQESSEAEEDMSVPRGPIAQWADGAISP  
NGHRPQLRAASGEDSTEEALAAAAAPWEGGPVPGVDTSPFAKSLGHSRGEADLFDSGDIFSTGTGSQSVE  
RTKPKAKIAENPANPPVGGKAKSPMFPALGEASSDDDLFQSAKPKPAKKTNPFLLEDEDDLFTDQKVKK  
NETKSSSQDVLTTQDIFEDDIFATEAIKPSQKTREKEKTLESNLFDDNIDIFADLTVKPKKEKSKKKVE  
AKSIFDDDDMDDIFSTGIQAKTTKPKSRSAAQAAPEPRFEHKVSNIFDDPLNAFGGQ


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	137
Concentration:	>0.05 µg/µL as determined by microplate BCA method



[View online »](#)

<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001162578</a>
<b>Locus ID:</b>	253725
<b>UniProt ID:</b>	<a href="#">Q9Y4E1</a> , <a href="#">B3KMC4</a>
<b>Cytogenetics:</b>	10q11.22
<b>RefSeq ORF:</b>	3735
<b>Synonyms:</b>	FAM21A; FAM21C; VPEF
<b>Summary:</b>	Acts at least in part as component of the WASH core complex whose assembly at the surface of endosomes inhibits WASH nucleation-promoting factor (NPF) activity in recruiting and activating the Arp2/3 complex to induce actin polymerization and is involved in the fission of tubules that serve as transport intermediates during endosome sorting. Mediates the recruitment of the WASH core complex to endosome membranes via binding to phospholipids and VPS35 of the retromer CSC. Mediates the recruitment of the F-actin-capping protein dimer to the WASH core complex probably promoting localized F-actin polymerization needed for vesicle scission (PubMed:19922874, PubMed:20498093, PubMed:22513087, PubMed:23331060). Via its C-terminus binds various phospholipids, most strongly phosphatidylinositol 4-phosphate (PtdIns-(4)P), phosphatidylinositol 5-phosphate (PtdIns-(5)P) and phosphatidylinositol 3,5-bisphosphate (PtdIns-(3,5)P2). Involved in the endosome-to-plasma membrane trafficking and recycling of SNX27-retromer-dependent cargo proteins, such as GLUT1 (PubMed:25278552). Required for the association of DNAJC13, ENTR1, ANKRD50 with retromer CSC subunit VPS35 (PubMed:24980502). Required for the endosomal recruitment of CCC and retriever complexes subunits COMMD1 and CCDC93 as well as the retrievere complex subunit VPS35L (PubMed:25355947, PubMed:28892079).[UniProtKB/Swiss-Prot Function]

**Product images:**

Coomassie blue staining of purified WASHC2C protein (Cat# [TP330681]). The protein was produced from HEK293T cells transfected with WASHC2C cDNA clone (Cat# [RC230681]) using MegaTran 2.0 (Cat# [TT210002]).