

Product datasheet for TP330683L

PCDH7 (NM_001173523) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protocadherin 7 (PCDH7), transcript variant d, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC230683 representing NM_001173523 Red=Cloning site Green=Tags(s)

MLRMRTAGWARGWCLGCCLLLPLSLSLAAAKQLLRYRLAEEGPADVIRGNVASDLGIVTGSGEVTFLES
GSEYLKIDNLTGELSTSERRIDREKLPQCQMIFDENECFLDFEVSIVIGPSQSWVDLFEGQVIVLDINDNT
PTFPSPVLTLTVEENRPVGTLYLLPTATDRDFGRNGIERYELLQEPGGGGSGGESRRAGAADSAPYPGGG
GNGASGGGSGSKRRLDASEGGGTNPGRSSVFELQVADTPDGEKQPQLIVKALDREQRDSYELTLRV
RDGGDPPRSSQAILRVLITDVNDNSPRFEKSVYEADLAENSAPGTPILQLRAADLDVGVNGQIEYVFGAA
TESVRRLLRLDETSWLSVLHRIDREEVNQLRFTVMARDRGQPPKTDKATVVLNIKDENDNVPSIEIRKI
GRIPLKDGVANVAEDVLVDTPIALVQVSDRDQGENGVVTCTVVGDPVFPQLKPASDTEGDQNKKKYFLHTS
TPLDYEATREFNVVIVAVDSGSPSLSSNNSLIVKVGDTNDNPPMFGQSVVEVYFENNIPGERVATLAT
DADSGKNAEIAYSLDSSVMGIFAIDPDSGDILVNTVLDREQTDRYEFKVNADKDGIPVLQGSTTVIVQVA
DKNDNDPKFMQDVFTFYKLENLQPNSPVGMVTVMDADKGRNAEMSLYIEENNNIFSIENDTGTIYSTMSF
DREHQTTYTFRVKAVDGGDPPRSATATVSLFVMDENDNAPTTLPKNISYLLPPSSNVRTVVATVLTAD
SDDGINADLNYSIVGGNPFKLFIDPTSGVSLVGLTKQKHYGLHRLVQVNDSGQPSQSTTLVHVFN
ESVSNATAIDSQIARSLHIPLTQDIAGDPSYEISKQRLSIVIGVVAGIMTVILILIVVMARYCRSKNKN
GYEAGKKDHEDFFTPQQHDKSKKPKKDKKNKSKQPLYSSIVTVEASKPNGQRYDSVNEKLSDSPSMGRY
RSVNGGPGSPDLARHYKSSSPLPTVQLHPQSPTAGKKHQAVQDLPPANTFVGAGDNISIGSDHCSEYSCQ
TNNKYSKQPFRRVTFVSVSQPDHPQGLQSCYDSGLEESETPSSKSSSGPRLGALPLPEDNYERTTDPG
SVGEAEHMENDSRPLPDVALTGKCTRECDEYGHSDSCWMPVRTSPERKKSQPKLSTFMPVDERGSQEKLA
NGEAAIMGDRNRNLLNKKLTSSYETFSAASFKNNEANPEDIPLTKTGEYKPSPVNTLTRREVYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	137
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001166994
Locus ID:	5099
UniProt ID:	O60245
Cytogenetics:	4p15.1
RefSeq ORF:	3765
Synonyms:	BH-Pcdh; BHPCDH; PPP1R120
Summary:	This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The gene encodes a protein with an extracellular domain containing 7 cadherin repeats. The gene product is an integral membrane protein that is thought to function in cell-cell recognition and adhesion. Alternative splicing yields isoforms with unique cytoplasmic tails. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified PCDH7 protein (Cat# [TP330683]). The protein was produced from HEK293T cells transfected with PCDH7 cDNA clone (Cat# [RC230683]) using MegaTran 2.0 (Cat# [TT210002]).