

Product datasheet for TP331069

MKRN2OS (NM_001195279) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC100129480 (LOC100129480), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC231069 representing NM_001195279 Red =Cloning site Green =Tags(s) MHCAEAGKALIKFNHCEKYIYSFSVPQCCPLCQQDLGSRKLEDAPVSIANPFTNGHQEKCSFLLRPTQGT FLREYDGRSDLHVGITNTNGVWYNYSAGVQRDGEWEESISIPLLQPNMYGMMEQWDKYLEDFSTSGAW LPHRYEDNHHNCYSYALTFINCVLMAEGRQQLDKGEFTEKYVVPRTLASKFITLYRAIREHGFYVTDGP QQQAQPPEGGLC TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	25.9
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001182208</u>
Locus ID:	100129480
UniProt ID:	<u>H3BPM6</u>



[View online »](#)

Cytogenetics: 3p25.2

RefSeq ORF: 669

Synonyms: C3orf83; MKRN2-AS1

Product images:



Coomassie blue staining of purified MKRN2OS protein (Cat# TP331069). The protein was produced from HEK293T cells transfected with MKRN2OS cDNA clone (Cat# [RC231069]) using MegaTran 2.0 (Cat# [TT210002]).