

Product datasheet for **TP331144**

IRF3 (NM_001197124) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Purified recombinant protein of Homo sapiens interferon regulatory factor 3 (IRF3), transcript variant 4, 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC231144 representing NM_001197124
Red=Cloning site **Green**=Tags(s)

MGTPKPRILPWLVSQLDLGQLEGVAWVNKSRTRFRIPWKHGLRQDAQQEDFGIFQAWAEATGAYVPRDK
PDLPTWKRNFERSALNRKEGLRLAEDRSKDPHDPHKIYEFVNSGVGDFSQPDTSPDTNGGGSTSDTQEDIL
DELLGNMVLAPLPDPGPPSLAVAPEPCQPLRSPSLDNPTFPNLGPESENPLKRLVLPGEDLITFTEGSG
RSPRYALWFCVGESWPQDQPWTKRLVMVKVPTCLRALVEMARVGGASSLENTVDLHISNSHPLSLTSDQ
YKAYLQDLVEGMDFQGPGES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 33.6
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_001184053](#)
Locus ID: 3661



[View online »](#)

UniProt ID: [Q14653](#), [E2GIM9](#), [E2GIM8](#)

Cytogenetics: 19q13.33

RefSeq ORF: 900


Synonyms: IIAE7

Summary: This gene encodes a member of the interferon regulatory transcription factor (IRF) family. The encoded protein is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other interferon-induced genes. The protein plays an important role in the innate immune response against DNA and RNA viruses. Mutations in this gene are associated with Encephalopathy, acute, infection-induced, herpes-specific, 7. [provided by RefSeq, Sep 2020]

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Cytosolic DNA-sensing pathway, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway

Product images:



Coomassie blue staining of purified IRF3 protein (Cat# TP331144). The protein was produced from HEK293T cells transfected with IRF3 cDNA clone (Cat# [RC231144]) using MegaTran 2.0 (Cat# [TT210002]).