

Product datasheet for TP331145L

C8orf58 (NM_001198827) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 8 open reading frame 58 (C8orf58), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC231145 representing NM_001198827 Red =Cloning site Green =Tags(s)
	MMGRRRAFAVDGRDGAGEGLARGCIVPGVTSTYRRIPDAAHGCSSWERGDKFRGVGREALFLKLASRDSG VEMAVGDSPLAALPGLSQDSLDFESSGSSEPPAQVGRLLASQKLGEVLERSRRLPTAPTSLSGQHRSLRL ASKPEREVPLGAGQQESMEADTDLEAGLEEEAVGGLPGAWACLPGQGLRYLEHLCLVLEQMARLQQLY L QLRIQRPPGDPGEEESTRAPLPSPLHTPGNRGQGPWELLSQTEHTGAKAASPPKVEVPSANPPRLPETPV EPTYHLPSSQGHKDRVQGPP
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.5
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001185756</u>



[View online »](#)


Locus ID: 541565

UniProt ID: [A0A087WX44](#)

Cytogenetics: 8p21.3

RefSeq ORF: 900

Product images:



Coomassie blue staining of purified C8orf58 protein (Cat# [TP331145]). The protein was produced from HEK293T cells transfected with C8orf58 cDNA clone (Cat# [RC231145]) using MegaTran 2.0 (Cat# [TT210002]).