

Product datasheet for **TP331187M**

C10orf81 (PLEKHS1) (NM_001193434) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens chromosome 10 open reading frame 81 (C10orf81), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC231187 representing NM_001193434 Red =Cloning site Green =Tags(s) MQSVQKMFKCHPDEVMSIRTTNREYFLIGHDREKIKDWVSFMSFRQDIKATQQNTEEELSLGNKRTLFI SSPLLGPSSTSEAVGSSSPRNLQDKHLMESQSPGFRQTHLQDLSEATQDVKEENHYLTPRSVLLLEDNI IASSDSGESIETDGPDQVSGRIECHYEPMESYFFKETSHEVSDSSKEEPQTLPETQDGDHLHQEQSGID WCLSPADVVAQTTNDQKGNIPDESQVEKLVFLSPPDVINYALALTEATGRICVSQWEGPPRLGCFCHGD HLLAVNDLKPQSLEEVSLFLTRSIQKEKLTIGRIPNSETFHAASCMCPSKCSAAPSQDKPRLNRP KRSPAIKKSQQKGARE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	41.4
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001180363</u>



[View online »](#)

Locus ID: 79949
UniProt ID: [Q5SXH7](#)
Cytogenetics: 10q25.3
RefSeq ORF: 1098
Synonyms: C10orf81; HEL185

Product images:



Coomassie blue staining of purified PLEKHS1 protein (Cat# [TP331187]). The protein was produced from HEK293T cells transfected with PLEKHS1 cDNA clone (Cat# [RC231187]) using MegaTran 2.0 (Cat# [TT210002]).