

Product datasheet for **TP331197L**

ANKRD63 (NM_001190479) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ankyrin repeat domain 63 (ANKRD63), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC231197 representing NM_001190479 Red =Cloning site Green =Tags(s) MLKPKDLCPRAGTRTFLEAMQAGKVHLARFVLDALDRSIIDCRAEQGRTPLMVAVGLPDPALRARFVRLLEQGAAVNLRDERGRTALSLACERGHLDVAVQLLVQFSGDPEAADSAGNSPVMWAAACGHGAVLEFLVRSFRRLLGLRLDRTNRAGLTALQLAAARGHGTCVQALTGPWGRAAAAAAARGSNSDSPGRPAPAASPEHRRPSRRLLPRLLARFARAAGGHGGEAGSAGKNSGRHRAQGSERPELGRSMLALGAVTEEEAARLRAGALMALPNSPQSSGTGRWRSQEVLEGAPPTLAQAPIGLSPHPEGPGSGRLGLRRRSTAPDIPSLVGEAPGPESGPELEANALSVSVPGNPWQAGTEAVVLAQR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	40.1
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001177408</u>
Locus ID:	100131244



[View online »](#)

UniProt ID: [C9JTQ0](#)
Cytogenetics: 15q15.1
RefSeq ORF: 1140

Product images:



Coomassie blue staining of purified ANKRD63 protein (Cat# [TP331197]). The protein was produced from HEK293T cells transfected with ANKRD63 cDNA clone (Cat# [RC231197]) using MegaTran 2.0 (Cat# [TT210002]).