

Product datasheet for TP331305M

SCHIP1 (IQJ-SCHIP1) (NM_001197114) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human IQJ-SCHIP1 readthrough (IQJ-SCHIP1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC231305 representing NM_001197114 Red =Cloning site Green =Tags(s)
	<p>MRLEELKRLQNPLEQVNDGKYSFENIQRAWREYLQRQEPLGKRSPSPSPSVSEKLSVVSMNTFSDSSTP DYREDGMDLGSDAGSSSSSRASSQSNSTKVTPCSECKSSSSPGGSLDLVSALEDYEEPPVYQKKVIDE WAPEEDGEEEEEDERDQRGYRDDRSPAREPGDVSARTRSGGGGGRSATTAMPPPVPNGNLHQHDPQDLR HNGNVVWVAGRPSCSRGPRRAIQKPQPAGGRRSGRGAAGGLCLQPPDGGTCVPEEPPVPPMDWEALEKHL AGLQFREQEVNRNQGQARTNSTSAQKNERESIRQKLALGFFDDGPGIYTSKSKGKPSLSSRLQSGMNLQ ICFVNDSGSDKSDADDSKTETSLDTPLSPMSKQSSSYSDRDTTEEESESLDDMDFLTRQKKLQAEAKMA LAMAQPMKMQVEVEKQNRKKSPVADLLPHMPHISECLMKRSLKPTDLRDMTIGQLQVIVNDLHSQIESL NEELVQLLLIRDELHTEQDAMLVDIEDLTRHAESQQKHMAEKMPAK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	59.6
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	NULL or Add: Recombinant proteins was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001184043
Locus ID:	100505385
UniProt ID:	Q9P0W5 , B3KU38
Cytogenetics:	3q25.32-q25.33
RefSeq ORF:	1608
Summary:	This locus represents naturally occurring read-through transcription from the neighboring IQ motif containing J (IQJ) and schwannomin interacting protein 1 (SCHIP1) genes. Alternative splicing results in multiple transcript variants that are composed of in-frame exons from each individual gene. The resulting fusion products are thought to be components of the multimolecular complexes of axon initial segments and nodes of Ranvier, and they may play a role in calcium-mediated responses. [provided by RefSeq, Oct 2010]

Product images:



Coomassie blue staining of purified IQJ-SCHIP1 protein (Cat# [TP331305]). The protein was produced from HEK293T cells transfected with IQJ-SCHIP1 cDNA clone (Cat# [RC231305]) using MegaTran 2.0 (Cat# [TT210002]).