

Product datasheet for **TP500064**

Gngt2 (NM_001038664) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 2 (Gngt2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200064 protein sequence Red=Cloning site Green=Tags(s)
	MAQDLSERELLRMEVEQLKKEVKNPRDLISKTGKEIKDYVEAQAGTDPLLKGIPEDKNPFKEKGTCLVLS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	7.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001033753
Locus ID:	14710
UniProt ID:	Q6101Z , Q6P8Y9 , A2A614
RefSeq Size:	586
Cytogenetics:	11 59.01 cM
RefSeq ORF:	210



[View online »](#)

Synonyms: AV096488

Summary: Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction. [UniProtKB/Swiss-Prot Function]