

## Product datasheet for **TP500076**

### Rpl38 (NM\_001048058) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ribosomal protein L38 (Rpl38), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200076 protein sequence Red=Cloning site Green=Tags(s)
	MPRKIEIEIKDFLLTARRKDAKSVKIKKNKDNVFKVRCSTRYLYTLVITDKEKAELKQSLPPGLAVKDLK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	8.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001041523</a>
Locus ID:	67671
UniProt ID:	<a href="#">Q9JII8</a> , <a href="#">Q52KP0</a>
RefSeq Size:	313
Cytogenetics:	11 E2
RefSeq ORF:	213



[View online »](#)

**Synonyms:** 0610025G13Rik; Rbt; Ts; Tss