

Product datasheet for **TP500092**

1110002E23Rik (BC021590) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA 1110002E23 gene (cDNA clone MGC:35768 IMAGE:5068609), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200092 protein sequence Red =Cloning site Green =Tags(s) MASSGVTVSAAGSASEASEVDPDNVGDWLRGVFRFATDRNDFRRNLILNLGLFAAGVWLARNLSDIDLMA P QPGV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	7.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	66119
UniProt ID:	Q9CQN3
RefSeq Size:	1059
Cytogenetics:	17 C



[View online »](#)

RefSeq ORF: 222

Synonyms: 1110002E23Rik; AI663970